



REF FYL-1S48-100

Rapid Response® Fentanyl Test Strip



Fentanyl is a synthetic opioid related to the phenylpiperidines that is commonly used for moderate to severe pain relief. Continued use of fentanyl can lead to opioid use disorder.

The Rapid Response® Fentanyl (FYL) Test Strip is an easy-to-use rapid lateral flow assay for the qualitative detection of fentanyl in at a cut-off of 200 ng/mL.

Kit Contents

- Individually packaged test strips
- Product insert


Product Details

- Built-in procedural control
- Results in 5 minutes
- Compact format for easy distribution

Detects Fentanyl Analogues including:

Carfentanil	Valeryl Fentanyl
Butyryl Fentanyl	Ocfentanil
p-Fluoro Fentanyl	3-Methyl Fentanyl
Acetyl Fentanyl	Remifentanil
Fentanyl	Sufentanil
Furanyl Fentanyl	

Testing Procedure



For illustrative purposes only. Refer to product insert for more details.

Results Interpretation

Rapid Response Results Interpretation Card
Fentanyl Test Strips

This card is not a replacement for the instructions for use. Please read the instructions prior to use and follow the directions carefully.

Positive Contains Fentanyl	
Negative Fentanyl could not be detected	
Invalid	

John Hopkins Research

Rapid Response® Fentanyl Test Strips showed the lowest detection limit and highest accuracy among other fentanyl detection methods when tested by John Hopkins Bloomberg School of Public Health.



VHC

The Vancouver Coastal Health Authority offers Rapid Response® Fentanyl Strips at their overdose prevention and supervised consumption sites and in their take-home drug testing kits.

Ordering Information

Product Code	Product Name	Contents
FYL-1S48-100	Fentanyl Single Drug Test Strip	100 Tests

To learn more contact your local **Sales Representative**, call us at **+1 888 339-9964**, or email us at **sales@btx.com**